



Recombinant Cupriavidus taiwanensis Thymidylate kinase (tmk)

Product Code	CSB-EP007226DZS
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B3R181
Product Type	Recombinant Protein
Immunogen Species	Cupriavidus taiwanensis (strain DSM 17343 / BCRC 17206 / CIP 107171 / LMG 19424 / R1) (Ralstonia taiwanensis (strain LMG 19424))
Purity	>85% (SDS-PAGE)
Sequence	MRGKFITFEG IDGAGKSTHI DWVADRLRAR ADIAGVVTTR EPGGTSLGED LRQILLHRKM HLETEALLMF AARREHIAEV IAPALERGKW VISDRFTDAT FAYQGGGRGL ATERLEVLEN WVQGS LQPD LLLFDVPLET ASARLAAART PDKFEAESRA FFQRTRDEYL RRAAQSPHRF RVIDATRSIA DIRDELEKIL ATI
Source	E.coli
Target Names	tmk
Protein Names	Recommended name: Thymidylate kinase EC= 2.7.4.9 Alternative name(s): dTMP kinase
Expression Region	1-203
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.