



Recombinant Human Cathepsin W (CTSW)

| | |
|--------------------------|---|
| Product Code | CSB-BP006205HU |
| Storage | Store at -20°C, for extended storage, conserve at -20°C or -80°C. |
| Uniprot No. | P56202 |
| Product Type | Recombinant Protein |
| Immunogen Species | Homo sapiens (Human) |
| Purity | ≥85% (SDS-PAGE) |
| Sequence | VPF SCDWRKVASA ISPIKDQKNC NCCWAMAAAG NIETLWRISF WDFVDVSVQE LLDCGRCDG CHGGFVWDAF ITVLNNSGLA SEKDYPFQGGK VRAHRCHPKK YQKVAWIQDF IMLQNNHRI AQYLATYGPI TVTINMKPLQ LYRKGVIKAT PTTCDPQLVD HSVLLVGFSG VKSEEGIWAE TVSSSQSQPQP PHPTPYWILK NSWGAQWGEK GYFRLHRGSN TCGITKFPPLT ARVQKPDMPK RVSCPP |
| Source | Baculovirus |
| Target Names | CTSW |
| Protein Names | Recommended name: Cathepsin W EC= 3.4.22.- Alternative name(s): Lymphopain |
| Expression Region | 128-376 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Full Length of Mature Protein |
| Target Details | This protein, a member of the peptidase C1 family, is a cysteine proteinase that may have a specific function in the mechanism or regulation of T-cell cytolytic activity. The encoded protein is found associated with the membrane inside the endoplasmic reticulum of natural killer and cytotoxic T-cells. Expression of this gene is up-regulated by interleukin-2. |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |