



# Recombinant *Drosophila grimshawi* Methylthioribose-1-phosphate isomerase (GH21593)

<b>Product Code</b>	CSB-YP014799DLL
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	B4JRX2
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Drosophila grimshawi</i> (Fruit fly) ( <i>Idiomya grimshawi</i> )
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MSLQSIKYKR GSLEILDQLL LPVVS KYLPV RGVEDGWKVI NKMQVRGAPA IAIVGCLSLA VEIYPEEFET KKSLRQEIEG KLN YLVSARP TAVNMKISAD ELITLANELS KNDDVTVANM KQRFLDATEA MLEKDIADNR AIGSNGAKAI LERVEATGS PGSASPVRVL THCNTGSLAT AGYGTALGVV RHLSELGKLE HVYCTETRPY NQGARLTAYE LVHEKLPATL VLDSMVAALF RVKNVAVVV GADRVAANGD TANKIGTYQI AVVAKHHGVP FYVAAPLTSI DLQIPSGEHI IIEVRPDREM THVGEHRIAA PGINCWNP AF DVTPASLITG IITEHGVFKP STLKDEIAKL IEL
<b>Source</b>	Yeast
<b>Target Names</b>	GH21593
<b>Protein Names</b>	Recommended name: Methylthioribose-1-phosphate isomerase Short name= M1Pi Short name= MTR-1-P isomerase EC= 5.3.1.23 Alternative name(s): S-methyl-5-thioribose-1-phosphate isomerase Translation initiation factor eIF-2B sub
<b>Expression Region</b>	1-363
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.