



Recombinant *Drosophila mojavensis* Flap endonuclease 1 (Fen1)

Product Code	CSB-MP008585DLY
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B4KNM1
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila mojavensis</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MGILGLSKLI ADLAPLAIRE SEIKNFFGRK VAIDASMCLY QFLIAVRSEG AQLATVNGDP TSHLMGMFYR TIRLLDNGIK PVYVFDGKPP DLKAGELAKR AERREEAEKA LKVATDAGDE AEIEKFNRRRL VRVTKEHSNE AKELLKLMGV PYVDAPCEAE AQCAALVKAG KVYATATEDM DALTFGSCKL LRYLTYSEAR KMPVKEFSYD KVLQGLELTS KEFIDLCILM GCDYCDSIKG IGPKRAIELI KTYRDIETIL ENIDTSKYIV PENWNYQRAR ELFVEPEVTD ASTIDLKWTA PDEDGLVQFL CGDRQFNEER VRNGARKLLK SKQSQTQVRL DSFFKALPSS PNATAAAKRK AEEIKKSANN KKAKTSGGSG AARGRRPK
Source	Mammalian cell
Target Names	Fen1
Protein Names	Recommended name: Flap endonuclease 1 Short name= FEN-1 EC= 3.1.-.- Alternative name(s): Flap structure-specific endonuclease 1
Expression Region	1-388
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.