



# Recombinant Porphyromonas gingivalis Lipoyl synthase (lipA)

<b>Product Code</b>	CSB-YP012927EXZ
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	B2RKU2
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Porphyromonas gingivalis (strain ATCC 33277 / DSM 20709 / CIP 103683 / JCM 12257 / NCTC 11834 / 2561)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MAQHVKKPEW LKIRLGGNEK FTETKSIVEG HCLHTICTSG KCPNMGECWS RGTATFMIGG DICTRACRFC NLTGPRKPL NEAEPHVAL SIKLMGLNHA VVTSDRDDL PDYGATHWVK TIQEIRRINS GVTLEVLIPD FKGRMDLVDM IIEASPDVIS HNLETVRRLT PSVRSVATYD TSLAVLRHIA QSGKMPAKTG MMLGLGETEE EILELMDDAL AAGVSVITIG QYLQPSRKNL PVVEYITPEQ FEHLRLVGIE KGFRTIESAP LVRSSYHAER HL
<b>Source</b>	Yeast
<b>Target Names</b>	lipA
<b>Protein Names</b>	Recommended name: Lipoyl synthase EC= 2.8.1.8 Alternative name(s): Lip-syn Short name= LS Lipoate synthase Lipoic acid synthase Sulfur insertion protein LipA
<b>Expression Region</b>	1-282
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.