



Recombinant Mouse Transcription factor E2F2 (E2f2)

Product Code	CSB-MP007341MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P56931
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	MLRAPRTLAP ATAQPTKSLP ALNPTELWPS GLSSPQLCPA TTATYYTSL YTQTVPSSVA LGTCLDAPH GPEGQIVRCA PAGRLPAKRK LDLEGIGRPT VPEFRTPKGK CIRVDGLPSP KTPKSPGEKT RYDTSGLLLT KKFIYLLSES EDGVLDLNWA AEVLVDVQKRR IYDITNVLEG IQLIRKSKN NIQWVGRELF EDPTRPSRQQ QLGQELKELM NAEQTLDQLI QSCSLSFKHL TEDNANKKLA YVTYQDIRAV GNFKEQTVIA VKAPPQTRLE VPDRAEENLQ IYLKSTQGPI EVYLCPEEGQ EPDSPAKEAL PSTSALSPIP DCAQPGCSTD SGIAETIEPS VLIPQPIPPP PPPPLPPAPS LVPLEATDNM LELSHPLLQQ TEDQFLSPIL AANSPLISFS PPLDQDEYLW GMDEGEGISD LFDSYDLGDL LIN
Source	Mammalian cell
Target Names	E2f2
Protein Names	Recommended name: Transcription factor E2F2 Short name= E2F-2
Expression Region	1-443
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Target Details	This protein is a member of the E2F family of transcription factors. The E2F family plays a crucial role in the control of cell cycle and action of tumor suppressor proteins and is also a target of the transforming proteins of small DNA tumor viruses. The E2F proteins contain several evolutionally conserved domains found in most members of the family. These domains include a DNA binding domain, a dimerization domain which determines interaction with the differentiation regulated transcription factor proteins (DP), a transactivation domain enriched in acidic amino acids, and a tumor suppressor protein association domain which is embedded within the transactivation domain. This protein and another 2 members, E2F1 and E2F3, have an additional cyclin binding domain. This protein binds specifically to retinoblastoma protein pRB in a cell-cycle dependent manner, and it exhibits overall 46% amino acid identity to E2F1.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the



contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.