



# Recombinant *Drosophila ananassae* Molybdopterin synthase catalytic subunit (GF17098)

<b>Product Code</b>	CSB-MP014706DKY
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	B3M268
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Drosophila ananassae</i> (Fruit fly)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MDHVKLVGDP IDISHIELL LDDGCGASSV FVGTTTRDNFQ GKKVVSLAYE AYDNMAIKEM NKICSDLRSK WPGLKHIVIH HRLGTVPVRE ASVIAASSP HRSEALESVS FAIDQLKTRV PIWKKEIYEG DHEAEWKENK ESIRAKQSLS SFNYSACPCQ VEESHEVPQK LVQIRVND AE LSKRLECFVN RKRDEINSQN IIDFKHPPAD SVPELSDSCA RTQSTVIKQE KSNCHIKVRR VNNRSFPQQM ELRPNYEEL NKLMSGRDGQ NEPSVKMRKS LYDSRLLTIE SYMGLAPKSE DDIFSRITL EDRILLLESI SPEYKHFTNI EPSSLEPPAP KKTRKKSYSY HELGMLIEKI KNE
<b>Source</b>	Mammalian cell
<b>Target Names</b>	Mocs2
<b>Protein Names</b>	Recommended name: Molybdopterin synthase catalytic subunit EC= 2.8.1.12 Alternative name(s): Molybdenum cofactor synthesis protein 2 large subunit Molybdenum cofactor synthesis protein 2B Short name= MOCS2B
<b>Expression Region</b>	1-363
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.