



Recombinant Human Ras-related protein Rab-25 (RAB25)

Product Code	CSB-YP019175HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P57735
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MGNGTEEDYN FVFKVVLIGE SGVGKTNLLS RFTRNEFSHD SRTTIGVEFS TRTVMLGTAA VKAQIWDTAG LERYRAITSA YYRGAVGALL VFDLTKHQTY AVVERWLKEL YDHAEATIVV MLVGNKSDLS QAREVPTTEA RMFAENGLL FLETSALDST NVELAFETVL KEIFAKVSKQ RQNSIRTNAI TLGSAQAGQE PGPGEKRACC
Source	Yeast
Target Names	RAB25
Protein Names	Recommended name: Ras-related protein Rab-25 Alternative name(s): CATX-8
Expression Region	1-210
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.