



Recombinant Human Sorting nexin-16 (SNX16)

| | |
|--------------------------|--|
| Product Code | CSB-BP022360HU |
| Storage | Store at -20°C, for extended storage, conserve at -20°C or -80°C. |
| Uniprot No. | P57768 |
| Product Type | Recombinant Protein |
| Immunogen Species | Homo sapiens (Human) |
| Purity | >85% (SDS-PAGE) |
| Sequence | <p>MATPYVPVPMPIGNSASSFTTNRNQRSSSFGSVSTSSNSSKGQLEDSNMGNF KQTSVPDQ MDNTSSVCSSPLIRTKFTGTASSIEYSTRPRDTEEQNPETVNWEDRPSTPTILG YEVMEE RAKFTVYKILVKKTPPEESWVVFRRYTDfSRLNDKLEMFPGFRLALPPKRWFK DNYNADF LEDRLGLQAFLQNLVAHKDIANCLAVREFLCLDDPPGPFDSLEESRAFCETLE ETNYRL QKELLEKQKEMESLKKLLSEKQLHIDTLENRIRTLSEPEESLDVSETEGEQILK VESSA LEVDQDVLDEESRADNKPCLSFSEPENAVSEIEVAEVAYDAEED</p> |
| Source | Baculovirus |
| Target Names | SNX16 |
| Protein Names | Recommended name: Sorting nexin-16 |
| Expression Region | 1-344 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Full length protein |
| Target Details | This gene encodes a member of the sorting nexin family. Members of this family contain a phox (PX) domain, which is a phosphoinositide binding domain, and are involved in intracellular trafficking. The function of this protein has not been determined. This gene results in three transcript variants encoding two distinct isoforms. |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |