



Recombinant Mouse Left-right determination factor 2 (Lefty2)

Product Code	CSB-YP012858MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P57785
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	<p>FSQ NFREVAGRFL MSETSTHLLV FGMEQRLPPN SELVQAVLRL FQEPVPR TAL RRFERLSPHS ARARVTIEWL RVREDGSNRT ALIDSRLVSI HESGWKAFDV TEAVNFWQQL SRPRQPLLLQ VSVQREHLGP GTWSAHKLV R FAAQGTPDGK GQGEPQLELH TLDLKDYGAQ GNCDPEVPVT EGTRCCRQEM YLDLQGMKWA ENWILEPPGF LTYECVGSC L QLPESLTIGW PFLGPRQCVA SEMTSLPMIV SVKEGGRTRP QVVSLPNMRV QTCSCASDGA LIPRGIDL</p>
Source	Yeast
Target Names	Lefty2
Protein Names	Recommended name: Left-right determination factor 2 Alternative name(s): Left-right determination factor B Protein lefty-2 Protein lefty-B
Expression Region	78-368
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	<p>This gene encodes a member of the TGF-beta family of proteins. The encoded protein is secreted and plays a role in left-right asymmetry determination of organ systems during development. The protein may also play a role in endometrial bleeding. Mutations in this gene have been associated with left-right axis malformations, particularly in the heart and lungs. Some types of infertility have been associated with dysregulated expression of this gene in the endometrium. Alternative processing of this protein can yield three different products. This gene is closely linked to both a related family member and a related pseudogene.</p>
Reconstitution	<p>We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.</p>



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.