



Recombinant Mouse T-cell leukemia/lymphoma protein 1A (Tcl1a)

Product Code	CSB-BP023313MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P56280
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	MATQRAHRAETPAHPNRLWIWEKHVYLDEFRRSWLPVVIKSNEKFQVILRQED VTLGEAM SPSQLVPYELPLMWQLYPKDRYRSCDSMYWQILYHIKFRDVEDMLLELIDSES NDE
Source	Baculovirus
Target Names	Tcl1a
Protein Names	Recommended name: T-cell leukemia/lymphoma protein 1A Alternative name(s): Oncogene TCL-1 Short name= Oncogene TCL1 Protein p14 TCL1
Expression Region	1-116
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.