



Recombinant *Drosophila ananassae* Phosphotriesterase-related protein (GF17207)

Product Code	CSB-MP018966DKY
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B3M070
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila ananassae</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MSTVQTVLGS ITPNLLGRTL THEHVALDFE HFYRPPPPDF ESELKAKISM STLGYVRLYP YSSKENVRFY DEEALEAAKK DILLYKKKHGG GSIVENSSYG LKRNLIIIVE LAKSTGVHFI AGTGHYIHAM QDARNASLTV EQLSDLYSKD IITGLEVNGK MVKCGFIGEV ASVYPIHDFE KHCICAAGEI QEVLGCGVSM HPRVTKAPF EIMRLYLEAG GRADKCVMSH LDRTIFDVDE LLEFAKLGCY IQYDLFGTEC SFYQLNTSVD MISDGQRIDN LMKLINEGLV DKLLMSHDIH TKHRLTSYGG HGYHHIYTNI LPRMFARGVT VEQVEQMTVT NPANWLAFNP
Source	Mammalian cell
Target Names	GF17207
Protein Names	Recommended name: Phosphotriesterase-related protein EC= 3.1.-.- Alternative name(s): Parathion hydrolase-related protein
Expression Region	1-350
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.