



Recombinant *Drosophila ananassae* KRR1 small subunit processome component homolog (dbe)

Product Code	CSB-MP012502DKY
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B3MM49
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila ananassae</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MSDSETEETK ISTEPVENAW AMKIPSFKPE DNPHEGMVEES SFATLFPKYR ERYLKEVWPL VEQCVAEHL KAELDLVEGS MIVKTSRKTW DPYIIKARD MIKLMARSVP FEQAKRVLQD DIGCDIIKIG NLVHKKEKFV KRRQRLIGPN GATLKSIELL TDCYVLVQGN TVAALGPYKG LQQVRDIVVE TMNNVHPIYN IKALMIKREL MKDPRLANED WSRFLPKFKN KNISKRKQPK VKKAKKEYTP FPPAQPEKSV GKQLASGEYF LNKEQKQAKR QQERSAKQAD AAKRQDERRN KDFVPPTEEA SSSRKRSSAE NASSKVDVKA LKAKLLKANK KKV
Source	Mammalian cell
Target Names	dbe
Protein Names	Recommended name: KRR1 small subunit processome component homolog Alternative name(s): KRR-R motif-containing protein 1 Protein dribble
Expression Region	1-343
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.