



Recombinant *Drosophila sechellia* Eukaryotic translation initiation factor 3 subunit K (GM15855)

Product Code	CSB-EP007540DMG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B4I7U3
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila sechellia</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MSHLVKMENG QSQTIQEMLG CIERYNPDHL KTLESYVQDQ AKNNTYDLEA NLAVLKLYQF NPHMLNFDIT YTILLKSLTS LPHTDFVMAK CLLLPPQQMKD ENVQTIIDLA DILERADFTL FWQRAEVNRN MFRHIAGFHD SIRKFSVSHVV GTTFQTIRKD LLKELGGIE DSTLESWIKR NGWKNQGGQL VIVAMQDDKI KTKNITEKIE FDNVGALMAQ CL
Source	<i>E. coli</i>
Target Names	GM15855
Protein Names	Recommended name: Eukaryotic translation initiation factor 3 subunit K Short name= eIF3k Alternative name(s): eIF-3 p25
Expression Region	1-222
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.