



Recombinant *Drosophila virilis* Flap endonuclease 1 (Fen1)

Product Code	CSB-MP008585DMP
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B4LM90
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila virilis</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MGILGLSKLI ADLAPQAIRE SEIKNFFGRK VAIDASMCLY QFLIAVRSEG AQLATVNGDP TSHLMGMFYR TIRLLDNGIK PVYVFDGQPP DLKSGELAKR AERREEAEKA LKAATEAGDE AEIDKFNRRR VRVTKEHSRE AKELLKLMGV PYVDAPCEAE AQCAALVKAG KVIATATEDM DALTFGSCCKL LRYLTYSEAR KMPVKEFSYD KVLQGLELETS KEFIDLCILM GCDYCESIKG VGPKRAIELI KSYRDIETIL ENIDTSKYPV PENWNYKRAR ELFIEPDVAD ASAILDKWTE PDEDGLVQFL CGDRQFNEER VRNGAKKLLK SKQAQTQVRL DSFFTTLPST PNASTSAKR KAEAEKKSAN NKKAKTSSAG RGRRPK
Source	Mammalian cell
Target Names	Fen1
Protein Names	Recommended name: Flap endonuclease 1 Short name= FEN-1 EC= 3.1.-.- Alternative name(s): Flap structure-specific endonuclease 1
Expression Region	1-386
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.