



Recombinant Human Putative RNA-binding protein 11 (RBM11)

Product Code	CSB-MP019401HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P57052
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MFPAQEEADR TVFVGNLEAR VREEILYELF LQAGPLTKVT ICKDREGKPK SFGFVCFKHP ESVSYAIAL NGIRLYGRPI NVQYRFGSSR SSEPANQSFE SCVKINSHNY RNEEMLVGRS SFPMQYFPIN NTSLPQEYFL FQKMQWHVYN PVLQLPYYEM TAPLPNSASV SSSLNHVPDL EAGPSSYKWT HQQPSDSDLY QMTAPLPNSA SVSSSLNHVP DLEAGPSSYK WTHQQPSDSD LYQMNKRKRQ KQTSDSDSST DNNRGNECSQ KFRKSKKKKR Y
Source	Mammalian cell
Target Names	RBM11
Protein Names	Recommended name: Putative RNA-binding protein 11 Alternative name(s): RNA-binding motif protein 11
Expression Region	1-281
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.