



Recombinant Human RWD domain-containing protein 2B (RWDD2B)

Product Code	CSB-MP020605HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P57060
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MKIELSMQPW NPGYSSEGAT AQETYTCPKM IEMEQAEAQL AELDLLASMF PGENELIVND QLAVAELKDC IEKKTMEGRS SKVYFTINMN LDVSDEKMAM FSLACILPFK YPAVLPEITV RSVLLSRSQQ TQLNTDLTAF LQKHCHGDVC ILNATEWVRE HASGYVSRDT SSSPTTGSTV QSVDLIFTRL WIYSHHIYNK CKRKNILEWA KELSLSGFSM PGKPGVVCVE GPQSACEEFW SRLRKLNWKR ILIRHREDIP FDGTNDETER QRKFSIFEK VFSVNGARGN HMDFGQLYQF LNTKGCGDVF QMFFGVEGQ
Source	Mammalian cell
Target Names	RWDD2B
Protein Names	Recommended name: RWD domain-containing protein 2B
Expression Region	1-319
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.