



Recombinant *Drosophila ananassae* Eukaryotic translation initiation factor 3 subunit H (eIF-3p40)

Product Code	CSB-BP007537DKY
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B3MP94
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila ananassae</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MANRGGRHAR NEESDNTINY VQCDGLAVMK MVKHCHEESS NMDLAQGALL GLVVDKCLEI TNCFFPKSG DETMDEEMYQ LTVMRRLRRV NVDHLHVGWY QSSDVGNLS LALLESQYHY QTSIEESVVV VYDTQKSSRG FLCLKAYRLT PQAIQMYKDG DFTPEAFRTL KVGYESLFAE IPIVIKNSPL TNIMMSELNE LLPEDKGHNF LDLGTASVLE NHMRSLIERV DELYQEAVRY NKYQQVVKQ DTEKHRALAK LAAENAVRTS KGEPTVPEEE VIKQFRPMPV PARLTATITS GQINTHAQHI AQFCSQSLAK LFITESLQNA KEAKETK
Source	Baculovirus
Target Names	eIF-3p40
Protein Names	Recommended name: Eukaryotic translation initiation factor 3 subunit H Short name= eIF3h Alternative name(s): Eukaryotic translation initiation factor 3 subunit 3
Expression Region	1-337
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.