



Recombinant *Drosophila erecta* DDRGK domain-containing protein 1 (GG24454)

Product Code	CSB-EP006597DLK-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B3P045
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila erecta</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	AP ETKPAAAPQR GVPQRAQEGV PRRAQIARNQ RNRLRQNAPA APAGQVAPAA GAPAARGDSD HEDEGQVDAD DARVPQGAVL DEKMGAKKRA KMEAKEQKRL QREQELHDRE QRKVKEAKEE AERKQQEDLE AEVERKRVEA ERLAKEERER KEHEEYLKMK AAFSVEEEGF EEGDADEQDS LLAGFIQYIR DNKVVVLEDL AVAFKLKTQQ VIDRIQELQA DGTLTGVIDD RGKFIYVSEE ELSAVAKFIK QRGRVSITEL AESSNNLINL TPISAGGEEA SS
Source	E.coli
Target Names	GG24454
Protein Names	Recommended name: DDRGK domain-containing protein 1
Expression Region	29-312
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.