



Recombinant *Cupriavidus taiwanensis* Lipoyl synthase (lipA)

Product Code	CSB-EP012927DZS-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B2AG37
Product Type	Recombinant Protein
Immunogen Species	<i>Cupriavidus taiwanensis</i> (strain DSM 17343 / BCRC 17206 / CIP 107171 / LMG 19424 / R1) (<i>Ralstonia taiwanensis</i> (strain LMG 19424))
Purity	>85% (SDS-PAGE)
Sequence	MSDALIAPNA SSSEAPQSPA EHYDPTRKQK SADKTARIPI KIVPAEKLKK PDWIRVKAAT GNSRFYEIKD ILRANNLVTV CEEASCPNIG ECFGKGTATF MIMGDKCTRR CPFCDVGHGR PDPLDVNEPG NLARTIAQLK LNYVVITSVD RDDLRDGGAQ HYVDCISQTR ELSPATRIEV LVPDFRGRDL KALDILQACP PDVMNHNMET VPRLYKQARP GADYAHSLKL LQEFKRRNPN VPTKSGLMVG LGETDEEILE VMRDMRAHDI DMLTIGQYLA PSNHHPVLR YVHPDTFKMF EEEAYKMGFT HAAVGAMVRS SYHADQQAHQ AGFA
Source	E.coli
Target Names	lipA
Protein Names	Recommended name: Lipoyl synthase EC= 2.8.1.8 Alternative name(s): Lip-syn Short name= LS Lipoate synthase Lipoic acid synthase Sulfur insertion protein LipA
Expression Region	1-334
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.