



Recombinant *Saccharomyces cerevisiae* COP9 signalosome complex subunit 5 (RRI1)

Product Code	CSB-EP005792SVP-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B3LH96
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain RM11-1a) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MSLSNKTVKE LRQLLKERYT VEDELTESIA LSSMRFKPSQ EPEFHALSQS SLLKTKLKQQ SSTDIPSYTH VLISKLSCEK ITHYAVRGGN IEIMGILMGF TLKDNIVVMD CFNLPVVGTE TRVNAQLESY EYMQYIDEM YNHNDGGDGR DYKGAKLNVV GWFHSHPGYD CWLSNIDIQT QDLNQRFDQP YVAIVVDPLK SLEDKILRMG AFRTIESKSD DNSATSYYEL ETIIFDSELN RALFETKLNL HCVIEDDESE QISLNLRLIDS MKQYSYLMDS KNRVTRIKLA TTSERVSNEN KKNIDYQNRS TRSQFCLNTQ RGDSTETSSF GSMFSGDNTS DVDMEDRNLT QFDSTDTSLC INGEPSIHVN RVERSSRSTD NFHNSKKRMN SNQEKCHDEG NDMLQRNVLE TDYARAKNRI LASKIKQYER LRFYKDTFTL
Source	E.coli
Target Names	RRI1
Protein Names	Recommended name: COP9 signalosome complex subunit 5 EC= 3.4.-.-
Expression Region	1-440
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.