



# Recombinant *Saccharomyces cerevisiae* Methionine aminopeptidase 2 (MAP2)

<b>Product Code</b>	CSB-BP013716SVP
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	B3LNM2
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Saccharomyces cerevisiae</i> (strain RM11-1a) (Baker's yeast)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	<p>MTDAEIENSP ASDLKELNLE NEGVEQQDQA KADESDPVES KKKKNKKKKK  KKSINVKKIEL LFPDGKYPEG AWMDYHQDFN LQRTTDEESR YLKRDLERAE  HWNDVRKGAE IHRRVRAIK DRIVPGMKLM DIADMIENTT RKYTGAENLL  AMEDPKSQGI GFPTGLSLNH CAAHFTPANAG DKTVLKYEDV MKVDYGVQVN  GNIIDSAFTV SFDPQYDNLL AAVKDATYTG IKEAGIDVRL TDIGEAIQEV  MESYEVEING ETYQVKPCRN LCGHSIAPYR IHGGKSVPIV KNGDTTKMEE  GEHFAIETFG STGRGYVTAG GEVSHYARSA EDHQVMPTLD SAKNLLKTID  RNFGTLPFGR RYLDRLGQEK YLFALNNLVR HGLVQDYPPPL NDIPGSYTAQ  FEHTILLHAH KKEVVSKGDD Y</p>
<b>Source</b>	Baculovirus
<b>Target Names</b>	MAP2
<b>Protein Names</b>	Recommended name: Methionine aminopeptidase 2 Short name= MetAP 2 EC= 3.4.11.18 Alternative name(s): Peptidase M 2
<b>Expression Region</b>	1-421
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.