



# Recombinant *Drosophila ananassae* Cytoplasmic tRNA 2-thiolation protein 2 (GF15629)

<b>Product Code</b>	CSB-EP003313DKY
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	B3MN57
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Drosophila ananassae</i> (Fruit fly)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MCSIGEDDFG DEGGTHAMVA GSVPSGIDLS PGEC SKCDIS SDELFQLNFR TPECRDCFLA YVRHKFRAAL GAAKVLPRDA EVLLVVDGSA ESLVLLDMLY YAQTQNTFKR LHCSARMVYL EDQPVQGRDP ADLESLEALS KRYEAFEFYV ITLGLPVGSL KLLKDYS PVA QESNELIQKL SQLRSLTSRQ DY LQQQRK NL IASVAQKLHC SHVFNSSISV DLATQLLTSI ALGRGGS AAL DVALLDDRLA AGIKLLRPLK DLTEQEVQFY VHAQRLQPLL MASSRYGQEL GDAASLQNL T SAFVANLQKN FSSTVSTVFR TGGKISESTH PEQAKCNHCQ SALDVELSDT LLAIEYSRSV SEAGVRLNEA IEDPEIISK R RMKLNDELCH ACSSIEADFK NGS
<b>Source</b>	E.coli
<b>Target Names</b>	GF15629
<b>Protein Names</b>	Recommended name: Cytoplasmic tRNA 2-thiolation protein 2
<b>Expression Region</b>	1-403
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.