



Recombinant Mouse G/T mismatch-specific thymine DNA glycosylase (Tdg)

| | |
|--------------------------|--|
| Product Code | CSB-EP023342MO |
| Storage | Store at -20°C, for extended storage, conserve at -20°C or -80°C. |
| Uniprot No. | P56581 |
| Product Type | Recombinant Protein |
| Immunogen Species | Mus musculus (Mouse) |
| Purity | ≥85% (SDS-PAGE) |
| Sequence | MDAEAARSYS LEQVQALYSF PFQQMMAEVP NMAVTTGQQV PAVAPNMATV TEQQVPEDAP VQEPAPEAPK RKRKRPRAAE PQEPVEPKKP ATSKKSGKST KSKEKQEKIT DAFKVKRKVD RFNGVSEAEL LTKTLPDILT FNLDIVIIGI NPGLMAAYKG HHYPGPGNHF WKCLFMSGLS EVQLNHMDDH TLPKYGIGF TNMVERTTPG SKDLSSKEFR EGGRILVQKL QKYQPRIAVF NGKCIYEIFS KEVFGVKVKN LEFGLQPHKI PDTETLCYVM PSSSARCAQF PRAQDKVHYY IKLKDLRDQL KGIERNADVQ EVQYTFDLQL AQEDAKKMAV KEEKYDPGYE AAYGGAYGEN PCNGEPCGIA SNGLTAHSAE PRGEAAPSDV PNGQWMAQSF AEQIPSFNNC GTREQEEESH A |
| Source | E.coli |
| Target Names | Tdg |
| Protein Names | Recommended name: G/T mismatch-specific thymine DNA glycosylase EC=3.2.2.29 Alternative name(s): C-JUN leucine zipper interactive protein JZA-3 Thymine-DNA glycosylase |
| Expression Region | 1-421 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Full length protein |
| Target Details | This protein belongs to the TDG/mug DNA glycosylase family. Thymine-DNA glycosylase (TDG) removes thymine moieties from G/T mismatches by hydrolyzing the carbon-nitrogen bond between the sugar-phosphate backbone of DNA and the mispaired thymine. With lower activity, this enzyme also removes thymine from C/T and T/T mispairings. TDG can also remove uracil and 5-bromouracil from mispairings with guanine. This enzyme plays a central role in cellular defense against genetic mutation caused by the spontaneous deamination of 5-methylcytosine and cytosine. This gene may have a pseudogene in the p arm of chromosome 12. |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final |



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.