



Recombinant *Drosophila sechellia* Eukaryotic translation initiation factor 3 subunit J (Adam)

Product Code	CSB-BP007539DMG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B4IMQ8
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila sechellia</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MADDWESAAD SDVVIRPTAA ASVNKWEGED EDEDIKDSWE DEEEKKDEEK PTKTEAPAKP KPNKALKAKL EQQALREEEA EAERLANLSP AEKLAEKRLR QKIQEASDLK HAQEAFGVTS TCGGLDAFNP ETKEEFKEFG ATLSWKVAQF RESEHFPQFV EDLVRSLCVN LSAADIKKVK MNVEVLHSEK LKLEKANAKK PAGKGKGKVT LRTENDDIDG YQKYGNDFTE DYDDFM
Source	Baculovirus
Target Names	eIF3j
Protein Names	Recommended name: Eukaryotic translation initiation factor 3 subunit J Short name= eIF3j
Expression Region	1-236
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.