



# Recombinant *Drosophila yakuba* Eukaryotic translation initiation factor 3 subunit E (eIF3-S6)

<b>Product Code</b>	CSB-EP007534DMR-B
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	B4PK98
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Drosophila yakuba</i> (Fruit fly)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MANFDLTRIN CQFLDRHLTF PLLEFLCGKE IYNQQELLEY ILETVNKTNM IDYTM DTRKR LNLSQEMPEE LVQRKAEVLA TLKQLQNEVA PIMKATDILK NGESMKDSKT FVNALQKDYN FKVEHLESAY KLAKYLYECG NYQESTSYLY FCLIVMSPND KNYLNVLWVGK LAEILTINW NTALEDLTRL RDIYIDANFS TIQALQQRTW LIHWSVLVFF NHPKGRDLII EMFLYKPLYL NAIQTMCPHI MRYLATAVVI NRTRRNALKD LIKVIQQESY TYRDPITEFL ECLYVNFDFE GARLKLHECQ TVILNDFIV ACLNEFVEDA RLMIFETFCR IHQCITISML ADKLNMPNE AECWIVNLIR NARLNAKIDS KLGHVVMGTQ PLSPYQQQLVE KIDSLSMRSE HLAGLIERKS KQKQNQESAD SWKYY
<b>Source</b>	E.coli
<b>Target Names</b>	eIF3-S6
<b>Protein Names</b>	Recommended name: Eukaryotic translation initiation factor 3 subunit E Short name= eIF3e Alternative name(s): Eukaryotic translation initiation factor 3 subunit 6
<b>Expression Region</b>	1-435
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.