



Recombinant *Drosophila yakuba* Eukaryotic translation initiation factor 3 subunit E (eIF3-S6)

Product Code	CSB-MP007534DMR
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B4PK98
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila yakuba</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MANFDLTRIN CQFLDRHLTF PLLEFLCGKE IYNQQELLEY ILETVNKTNM IDYTM DTRKR LNLSQEMPEE LVQRKAEVLA TLKQLQNEVA PIMKATDILK NGESMKDSKT FVNALQKDYN FKVEHLESAY KLAKYLYECG NYQESTSYLY FCLIVMSPND KNYLNVLWVGK LAAEILTLNW NTALEDLTRL RDIYIDANFS TIQALQQRTW LIHWSVLVFF NHPKGRDLII EMFLYKPLYL NAIQTMCPHI MRYLATAVVI NRTRRNALKD LIKVIQQESY TYRDPITEFL ECLYVNFDFE GARLKLHECQ TVILNDFIV ACLNEFVEDA RLMIFETFCR IHQCITISML ADKLNMPNE AECWIVNLIR NARLNAKIDS KLGHVVMGTQ PLSPYQQQLVE KIDSLSMRSE HLAGLIERKS KQKQNQESAD SWKY Y
Source	Mammalian cell
Target Names	eIF3-S6
Protein Names	Recommended name: Eukaryotic translation initiation factor 3 subunit E Short name= eIF3e Alternative name(s): Eukaryotic translation initiation factor 3 subunit 6
Expression Region	1-435
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.