



# Recombinant Mouse Toll-like receptor 8 (Tlr8), partial

<b>Product Code</b>	CSB-YP023607MO
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P58682
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	IFHKANYSRSYPCDEIRHNSLVIAECNHRQLHEVPQTIGKYVTNIDLSDNAITHIT KESFQKLQNLTKIDLHNHNAKQQHPNENKNGMNITEGALLSLRNLTVLLLEDNQL YTIPAGLPESLKELSLIQNNIFQVTKNNTFGLRNLERLYLGWNCYFKCNQTFKV EDGAFKNLIHLKVL SLSFNLFYVPPKLPSSLRKLFSLNAKIMNITQEDFKLEN LTL DLSGNCPRCYNAPFPCTPCKENSSIIHPLAFQSLTQLLYLNL SSTSRLTI PSTWFENLSNLKELHLEFNLYVQEIASGAFLTKLPSLQILDLSFNFYKEYLQFI NISSNFSKLRSLKHLRGYVFRELK KKHFEHLQSLPNLATINLGINFIEKIDFKA FQNF SKLDVIYLSGNRIASVLDGTDYSSWRNRLRKPLSTDDDEFDPHVNFYHS TKPLIKPQCTAYGKALDLSLNNIFIIGKSQFEGFQDIACLNLSFNANTQVFNTE FSSMPHIKYLDLTNNRLDFDDNNAFSDLHDLEVL DLSHNAHYFSIAGVTHRLGF IQNLINLRVNL SHNGIYTLTEESELKSISL KELVFSGNRLDRLWNANDGKYWSI FKSLQNLIRLDLSYNNLQQIPNGAFLNLPQSLQELLISGNKLRFFNWTL LQYFPH LHLLDLSRNELYFLPNCLSKFAHSLETLLLSHNHF SHLPSGFLSEARNLVHLDLS FNTIKMINKSSLQTKMKTNLSILELHGNYFDCTCDISDFRSWLDENLNITIPKLVN VICSNPGDQKSKSIMSLDLTTCVSDTT
<b>Source</b>	Yeast
<b>Target Names</b>	Tlr8
<b>Protein Names</b>	Recommended name: Toll-like receptor 8 Alternative name(s): CD_antigen=CD288
<b>Expression Region</b>	24-818
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Partial
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.