



Recombinant Human Tubulin polymerization-promoting protein family member 2 (TPPP2)

Product Code	CSB-BP024116HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P59282
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MASEAEKTFH RFAAFGESSS SGTEMNKNF SKLCKDCGIM DGKTVTSTDV DIVFSKVKAK NARTITFQQF KEAVKELGQK RFKKGKSPDEV LENIYGLMEG KDPATTGATK ATTVGAVDRL TDTSKYTGTH KERFDESGKG KGIAGREEMT DNTGYVSGYK GSGTYDKKTK
Source	Baculovirus
Target Names	TPPP2
Protein Names	Recommended name: Tubulin polymerization-promoting protein family member 2 Alternative name(s): Protein p25-beta TPPP/p18
Expression Region	1-170
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.