



Recombinant Escherichia coli O157:H7 Peptide chain release factor 1 (prfA)

Product Code	CSB-YP015197EOE
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B5YXM9
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O157:H7 (strain EC4115 / EHEC)
Purity	>85% (SDS-PAGE)
Sequence	MKPSIVAKLE ALHERHEEVQ ALLGDAQTIA DQERFRALSR EYAQLSDVSR CFTDWQQVQE DIETAQMMLD DPEMREMAQD ELREAKEKSE QLEQQQLQVLL LPKDPDDERN AFLEVRAGTG GDEAALFAGD LFRMYSRYAE ARRWRVEIMS ASEGEHGGYK EIIAKISGDG VYGRCLKFESG GHRVQRVPAT ESQGRIHTSA CTVAVMPELP DAELPDINPA DLRIDTFRSS GAGGQHVNTT DSAIRITHLP TGIVVECQDE RSQHKNKAKA LSVLGARIHA AEMAKRQQAE ASTRRNLLGS GDRSDRNRTY NFPQGRVTDH RINLTLYRLD EVMEGKLDML IEPIIQEHQA DQLAALSEQE
Source	Yeast
Target Names	prfA
Protein Names	Recommended name: Peptide chain release factor 1 Short name= RF-1
Expression Region	1-360
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.