



Recombinant *Cryptosporidium muris* Flap endonuclease 1 (FEN1)

Product Code	CSB-YP008585DYT
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B6AFP1
Product Type	Recombinant Protein
Immunogen Species	<i>Cryptosporidium muris</i> (strain RN66)
Purity	>85% (SDS-PAGE)
Sequence	<p> MGIKGLARFL VDNAPKSIQQ QSIDSLIGRI IAIDASMWMY QFLAAIREGS QWGNLTNEAG ESTSHISGML SRTIRLLEAG IKPVFVFDGE PPELKMEELM KRKERREKAQ QELEKAQEEG DTETIRKQLI RTIKITKEQS DDVKYMLKLL GIPVIEATSE AEAQCAELCK EGLVYGVATE DADSLTFGTP LVIRHLNFSD GKFSDSKGNS KNSLQVIKLS TVLSDLGLSM QQFVDLCILC GCDYCGTIRG IGALTAYKLL KKHKDIETII NELDKNKHPL PMSFDYIKVR ELFTNPNVIK AKEFRDKLKW TNPNLEGLSE WLIKQQNFSE VRVANYCTRI KSKGKTAQT SLDSFFTLTP KKTSENKNKE KEITKVQEVK ETISNINKDT RKDSVSEWGE IKTKIETKE DKLDFERRNK IIQEDKLNSK KEFDYNEKLQ CNSTAENTPI IDNKAKVLSI KKRRINRVNY DE </p>
Source	Yeast
Target Names	FEN1
Protein Names	Recommended name: Flap endonuclease 1 Short name= FEN-1 EC= 3.1.-.- Alternative name(s): Flap structure-specific endonuclease 1
Expression Region	1-472
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.