



# Recombinant Human Transcription factor SOX-8 (SOX8)

<b>Product Code</b>	CSB-BP022436HU
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P57073
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	<p>MLDMSEARSQPPCSPSGTASSMSHVEDSDSDAPPSPAGSEGLGRAGVAVGG            ARGDPAEAA            DERFPACIRDAVSQVLKGYDWSLVPMPVRGGGGGALKAKPHVKRPMNAFMV            WAQAARRKL            ADQYPHLHNAELSKTLGKLWRLLESEKRPFVEEAERLRVQHKKDHPDYKYQ            PRRRKSAAK            AGHSDSDSGAELGPHPGGGAVYKAEAGLGDGHHHGDHTGQTHGPPTPPTTP            KTELQQAGA            KPELKLEGRRPVDSGRQNIDFSNVDISELSSEVMGTMDAFDVHEFDQYLPLGG            PAPPEPG            QAYGGAYFHAGASPVWAHKSAPSASASPTETGPPRPHIKTEQPSPGHYGDQ            PRGSPDYGS            CSGQSSATPAAPAGPFAGSQGDYGDLQASSYYGAYPGYAPGLYQYPCFHSP            RRPYASPLL NGLALPPAHSPTSHWDQPVYTTLTRP</p>
<b>Source</b>	Baculovirus
<b>Target Names</b>	SOX8
<b>Protein Names</b>	Recommended name: Transcription factor SOX-8
<b>Expression Region</b>	1-446
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full length protein
<b>Target Details</b>	This gene encodes a member of the SOX (SRY-related HMG-box) family of transcription factors involved in the regulation of embryonic development and in the determination of the cell fate. The encoded protein may act as a transcriptional activator after forming a protein complex with other proteins. This protein may be involved in brain development and function. Haploinsufficiency for this protein may contribute to the mental retardation found in haemoglobin H-related mental retardation (ART-16 syndrome).
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.