



Recombinant Mouse Matrix metalloproteinase-20 (Mmp20)

Product Code	CSB-YP014667MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P57748
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	YRLF PGEPKWKKNI LTYRISKYTP SMSPTVDKA IQMALHAWST AVPLNFVRIN SGEADIMISF ETGDHGDSYP FDGPRGTLAH AFAPGEGGLGG DTHFDNAEKW TMGTNGFNLF TVAAHEFGHA LGLGHSTDPS ALMYPTYKYQ NPYRFHLPKD DVKGIQALYG PRKIFPGKPT MPHIPP HKPS IPDLCDSSSS FDAVTMLGKE LLFFKDRIFW RRQVHLPTGI RPSTITSSFP QLMSNVDAAY EVAERGIAFF FKGP HYWVTR GFHMQGPRT IYDFGFPRHV QRIDAAVYLK EPQKTLFFVG EEYYSYDERK KKMEKDYPKN TEEEFSGVSG HIDA AVELNG YIYFFSGRKT FKYDTEKEDV VSVVKSSSWI GC
Source	Yeast
Target Names	Mmp20
Protein Names	Recommended name: Matrix metalloproteinase-20 Short name= MMP-20 EC= 3.4.24.- Alternative name(s): Enamel metalloproteinase Enamelysin
Expression Region	107-482
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	Proteins of the matrix metalloproteinase (MMP) family are involved in the breakdown of extracellular matrix in normal physiological processes, such as embryonic development, reproduction, and tissue remodeling, as well as in disease processes, such as arthritis and metastasis. Most MMP s are secreted as inactive proproteins which are activated when cleaved by extracellular proteinases. This protein degrades amelogenin, the major protein component of dental enamel matrix, and so the protein is thought to play a role in tooth enamel formation. A mutation in this gene, which alters the normal splice pattern and results in premature termination of the encoded protein, has been associated with amelogenesis imperfecta. This gene is part of a cluster of MMP genes that localizes to chromosome 11q22.3.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.