



Recombinant *Penicillium chrysogenum* Flap endonuclease 1 (fen1)

Product Code	CSB-EP008585ETH-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B6HEM2
Product Type	Recombinant Protein
Immunogen Species	<i>Penicillium rubens</i> (strain ATCC 28089 / DSM 1075 / NRRL 1951 / Wisconsin 54-1255) (<i>Penicillium chrysogenum</i>)
Purity	>85% (SDS-PAGE)
Sequence	MGIKHLYQVI SENAPDAIKT GEIKNHFGGRK VAIDASMSIY SFLIAVRSEG QQLMSESGET TSHLMGMFYR TLRMVDNGIK PLYVFDGAPP KLKSGELAKR VARKAEATEA HEEAKETGTA EDIEKFSRRT VRVTREHNAE CKKLELMGI PFINAPTEAE AQCAVLARAG KVYAAASEDM DTLCFESPIL LRHLTFSEQR KEPIQEIHL D RALEGLDMDR AQFIDLCILL GCDYLEPIPK VGATTALSLI KEHKSLEKVL EFMKNDPKKK FVVPEDWPYE DARELFTNPD VRPADHPECD FKWEAPNVEG LIEFLVGDKG FNEDRVRNGA ARLSKHLKTA QQRLEGFFK PVARTEDEKA SLKRKHDEKL QQQKKKKKED AKAKKEAKAK PRGAA
Source	<i>E. coli</i>
Target Names	fen1
Protein Names	Recommended name: Flap endonuclease 1 Short name= FEN-1 EC= 3.1.-.- Alternative name(s): Flap structure-specific endonuclease 1
Expression Region	1-395
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.