



Recombinant Acidithiobacillus ferrooxidans Lipoyl synthase (lipA)

Product Code	CSB-BP012927AWH
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B7J7N7
Product Type	Recombinant Protein
Immunogen Species	Acidithiobacillus ferrooxidans (strain ATCC 23270 / DSM 14882 / CIP 104768 / NCIMB 8455) (Ferrobacillus ferrooxidans (strain ATCC 23270))
Purity	>85% (SDS-PAGE)
Sequence	MDSQPQSKKA ARGADKTARN PIPIIPAPTE RLPKPQWLRV RSPLSPEVDQ LKKILRDAAL HTVCEEASCP NLGECFGGGT ATFMILGDIC TRRCPCFDVA HGRPEAPDPL ESVHLARTVA SMKLRV VIT SVDRDDL RDG GARHFAEVIS ALREHCP ELH VEILVPDFRG RVANALAAFR ET PPDVFNHN LETVPRLYLQ ARPGADYHHS LQLLA AFKAQ HPQIPTK SGL MLGLGEEIDE IRGVMRDLRQ AGCELLTIGQ YLAPSRHHLP VARFVPPEEF QQLQRDGMAM GFRHVASGPL VRSSYHAERS FHEIGGD
Source	Baculovirus
Target Names	lipA
Protein Names	Recommended name: Lipoyl synthase EC= 2.8.1.8 Alternative name(s): Lip-syn Short name= LS Lipoate synthase Lipoic acid synthase Sulfur insertion protein LipA
Expression Region	1-317
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.