



# Recombinant Mouse Cathepsin W (Ctsw)

<b>Product Code</b>	CSB-EP006205MO-B
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P56203
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	VPRTC DWRKAKNIIS SVKNQGSCKC CWAMAAADNI QALWRIKHQQ FVDVSVQELL DCERCGNGCN GGFVWDAYLT VLNNGLASE KDYPFQGDRK PHRCLAKKYK KVAWIQDFTM LSNNEQIAIH YLAVHGPIV TINMKLLQHY QKGVIKATPS SCDPRQVDHS VLLVGFGKEK EGMQTGTVLS HSRKRHSSP YWILKNSWGA HWGEKGYFRL YRGNNTCGVT KYPFTAQVDS PVKKARTSCP P
<b>Source</b>	E.coli
<b>Target Names</b>	Ctsw
<b>Protein Names</b>	Recommended name: Cathepsin W EC= 3.4.22.- Alternative name(s): Lymphopain
<b>Expression Region</b>	126-371
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Target Details</b>	This protein, a member of the peptidase C1 family, is a cysteine proteinase that may have a specific function in the mechanism or regulation of T-cell cytolytic activity. The encoded protein is found associated with the membrane inside the endoplasmic reticulum of natural killer and cytotoxic T-cells. Expression of this gene is up-regulated by interleukin-2.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.