



Recombinant Escherichia coli Galactokinase (galk)

Product Code	CSB-MP009198ENW
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B6I7R1
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain SE11)
Purity	>85% (SDS-PAGE)
Sequence	MSLKEKTQSL FANAFGYPAT HTIQAPGRVN LIGEHTDYND GFVLP CAIDY QTVISCAPRD DRKVRVMAAD YENQLDEFSL DAPIVAHENY QWANYVRGVV KHLQLRNNSF GGVDMVISGN VPQGAGLSSS ASLEVAVGTV LQQLYHLPLD GAQIALNGQE AENQFVGCNC GIMDQLISAL GKGDHALLID CRSLGKAVS MPKGVAVVII NSNFKRTL VG SEYNTRREQC ETGARFFQQP ALRDVTIEEF NAVAHELDPI VAKRVRHILT ENARTVEAAS ALEQGD LKRM GELMAESHAS MRDDFEITVP QIDTLVEIVK AVIGDKGGVR MTGGGFGGCI VALIPEELVP AVQQAVAEQY EAKTGIKETF YVCKPSQGAG QC
Source	Mammalian cell
Target Names	galk
Protein Names	Recommended name: Galactokinase EC= 2.7.1.6 Alternative name(s): Galactose kinase
Expression Region	1-382
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.