



# Recombinant Human Gasdermin-D (GSDMD)

<b>Product Code</b>	CSB-YP009956HU
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P57764
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MGSAFERVVR RVVQELDHGG EFIPVTSLQS STGFQPYCLV VRKPSSSWFW KPRYKCVNLS IKDILEPDAA EPDVQRGRSF HFYDAMDGQI QGSVELAAPG QAKIAGGAAV SDSSSTSMNV YLSVDPNTW QTLLHERHLR QPEHKVLQQL RSRGDNVYVV TEVLQTQKEV EVTRTHKREG SGRFSLPGAT CLQGEGQGHL SQKKTVTIPS GSTLAFRVAQ LVIDSDLVL LFPDKKQRTF QPPATGHKRS TSEGAWPQLP SGLSMMRCLH NFLTDGVP AE GAFTEDFQGL RAEVETISKE LELLDRELCQ LLEGLEGLV RDQLALRALE EALEQGQSLG PVEPLDGPAG AVLECLVLSS GMLVPELAIP VVYLLGALTM LSETQHKLLA EALESQTLLG PLELVGSLLE QSAPWQERST MSLPPGLLGN SWGEGAPAWV LLDECGLELG EDTPHVCWEP QAQGRMCALY ASLALLSGLS QEPH
<b>Source</b>	Yeast
<b>Target Names</b>	GSDMD
<b>Protein Names</b>	Recommended name: Gasdermin-D Alternative name(s): Gasdermin domain-containing protein 1
<b>Expression Region</b>	1-484
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full length protein
<b>Target Details</b>	Gasdermin D is a member of the gasdermin family. Members of this family appear to play a role in regulation of epithelial proliferation. Gasdermin D has been suggested to act as a tumor suppressor. Alternatively spliced transcript variants have been described.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.