



Recombinant *Taeniopygia guttata* Ragulator complex protein LAMTOR2 (LAMTOR2)

Product Code	CSB-EP020053TKW
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B5FYY5
Product Type	Recombinant Protein
Immunogen Species	<i>Taeniopygia guttata</i> (Zebra finch) (<i>Poephila guttata</i>)
Purity	≥85% (SDS-PAGE)
Sequence	MLRPKALTQV LSQANTGGVQ STLLLNNEGS LLAYSGYGDT DARVTAAIAS NIWVAYDKNG HQAFNEDNLK FILMDCMEGR VAITRVANLL LCMYAKETVG FGMLKAKAQA LVQYLEEPLT QVAAS
Source	E.coli
Target Names	LAMTOR2
Protein Names	Recommended name: Ragulator complex protein LAMTOR2 Alternative name(s): Late endosomal/lysosomal adaptor and MAPK and MTOR activator 2
Expression Region	1-125
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.