



Recombinant Escherichia coli O157:H7 NAD-dependent malic enzyme (maeA)

Product Code	CSB-YP013633EOE
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B5Z1T9
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O157:H7 (strain EC4115 / EHEC)
Purity	≥85% (SDS-PAGE)
Sequence	MEPKTKKQRS LYIPYAGPVL LEFPLLKGS AFSMEERRNF NLLGLLPEVV ETIEEQAERA WIQYQGFKTE IDKHIYLRNI QDNETLFYR LVNNHLDEMM PVIYTPTVGA ACERFSEIYR RSRGVFISYQ NRHNMDLILQ NVPNHNIKVI VVTDGERILG LGDQGIGGMG IPIGKLSLYT ACGGISPAYT LPVVLDVGTN NQQLLNPLY MGWRNPRITD DEYYEFVDEF IQAVKQRWPD VLLQFEDFAQ KNAMPLNRY RNEICSFNDD IQGTAAVTVG TLIAASRAAG GQLSEKKIVF LGAGSAGCGI AEMIIAQTR EGLSEEAARQ KVMVDRFGL LTKMPNLLP FQTKLVQKRE NLSDRDSD VLSLLDVVRN VKPDILIGVS GQTGLFTEEI IREMHKHCPR PIVMPLSNPT SRVEATPQDI IAWTEGNALV ATGSPFNPVV WKDKIYPIAQ CNNAFIFPGI GLGVIASGAS RITDEMLMSA SETLAQYSPL VLNGEGLVLP ELKDIQKVSRAIAFAVGKMA QQQGVAVKTS AEALQQAIDD NFWQAEYRDY RRTSI
Source	Yeast
Target Names	maeA
Protein Names	Recommended name: NAD-dependent malic enzyme Short name= NAD-ME EC= 1.1.1.38
Expression Region	1-565
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.