



Recombinant Zea mays Inosine triphosphate pyrophosphatase

Product Code	CSB-MP011907ZAX
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B6TNW8
Product Type	Recombinant Protein
Immunogen Species	Zea mays (Maize)
Purity	>85% (SDS-PAGE)
Sequence	MSAAARVLPK AVTFVTGNAK KLEEVRAILG SSVPFQSLKL DLPELQGEPE YISKEKARIA ASQVNGPVLV EDTCLCFNAL KGLPGPYIKW FLEKIGHEGL NLLKAYEDK SAFAMCIFSL ALPGGEEPIT FVGKTAGKIV PARGPNYFGW DPVFQPDGFE QTYAEMPKSV KNNISHRGKA LALVKEHFAS ASYTVQSNDST
Source	Mammalian cell
Protein Names	Recommended name: Inosine triphosphate pyrophosphatase Short name= ITPase Short name= Inosine triphosphatase EC= 3.6.1.19 Alternative name(s): Non-canonical purine NTP pyrophosphatase Non-standard purine NTP pyrophosphatas
Expression Region	1-201
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.