



Recombinant Mouse SET domain-containing protein 4 (Setd4)

Product Code	CSB-YP021105MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P58467
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	MQRRRGRTERARKRRRRSSGSRAVNESYRSEFIELRKWLKERKFEDTDLVPA SFPGTGRG LMSKASLQEGQVMISLPESCLLTDTVIRSSLGPYIKKWKPPVSPLLALCTFLVS EKHAG CRSLWKSYPDILPKSYTCPVCLEPEVVDLLPSPLKAKAEQRRARVQDLFTSAR GFFSTLQ PLFAEPVDSVFSYRAFLWAWCTVNTRAVYLRSSRQECLSAEPDTCALAPFLDL LNHSPHV QVKAAFNEKTRCYEIRTASRCRKHQEVFICYGPHDNQRLLEYGFSVRNPHA CVPVSAD MLVKFLPAADKQLHRKITILKDHGFTGNLTFGWDGPSWRLLTALKLLCLEAERF TSWKKV LLGEISDTNEKTSLGVAQKICSDVIEETHAVLRKVSDMKEGTVSLRNQLSLVE ALRMEE LRILQASAEILSGLLAPFS
Source	Yeast
Target Names	Setd4
Protein Names	Recommended name: SET domain-containing protein 4
Expression Region	1-439aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.