



Recombinant Escherichia fergusonii Galactokinase (galk)

| | |
|--------------------------|--|
| Product Code | CSB-MP009198EOR |
| Storage | Store at -20°C, for extended storage, conserve at -20°C or -80°C. |
| Uniprot No. | B7LK02 |
| Product Type | Recombinant Protein |
| Immunogen Species | Escherichia fergusonii (strain ATCC 35469 / DSM 13698 / CDC 0568-73) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MSLKEKTQSL FANAFGYPAT HTIQAPGRVN LIGEHTDYND GFVLPCAIDY QTVISCAPRD DRNVRVMAAD YENQLDEFSL DAPIVAHENY QWANYVRGVV KHLQLRNNSF GGVDMVISGN VPQGAGLSSS ASLEVAVGTV LQQLYHLPLD GAQIALNGQE AENQFVGCNC GIMDQLISAL GKKNHALLID CRSLGKAVS MPKGVAVVII NSNFKRTLVG SEYNTRREQC ETGARFFQQP ALRDVTIEEF NAVAHELDPI VAKRVRHILT ENARTVEAAS ALEQGD LKRM GELMAESHAS MRDDFEITVP QIDTLVEIVK AVIGDKGGVR MTGGGFGGCI VALIPEELVP AVQQAVAEQY EAKTGKETF YVCKPSQGAG QC |
| Source | Mammalian cell |
| Target Names | galk |
| Protein Names | Recommended name: Galactokinase EC= 2.7.1.6 Alternative name(s): Galactose kinase |
| Expression Region | 1-382 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |