



Recombinant Escherichia coli O8 Peptide chain release factor 1 (prfA)

| | |
|--------------------------|--|
| Product Code | CSB-BP015197EOO |
| Storage | Store at -20°C, for extended storage, conserve at -20°C or -80°C. |
| Uniprot No. | B7LXC6 |
| Product Type | Recombinant Protein |
| Immunogen Species | Escherichia coli O8 (strain IA11) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MKPSIVAKLE ALHERHEEVQ ALLGDAQTIA DQERFRALSR EYAQLSDVSR CFTDWQQVQE DIETAQMMLD DPEMREMAQD ELREAKEKSE QLEQQQLQVLL LPKDPDDERN AFLEVRAGTG GDEAALFAGD LFRMYSRYAE ARRWRVEIMS ASEGEHGGYK EIIAKISGDG VYGRCLKFESG GHRVQRVPAT ESQGRIHTSA CTVAVMPELP DAELPDINPA DLRIDTFRSS GAGGQHVNTT DSAIRITHLP TGIVVECQDE RSQHKNKAKA LSVLGARIHA AEMAKRQQAE ASTRRNLLGS GDRSDRNRTY NFPQGRVTDH RINLTLYRLD EVMEGKLDML IEPIIQEHQA DQLAALSEQE |
| Source | Baculovirus |
| Target Names | prfA |
| Protein Names | Recommended name: Peptide chain release factor 1 Short name= RF-1 |
| Expression Region | 1-360 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |