



Recombinant Bovine Lymphocyte antigen 96 (LY96)

| | |
|--------------------------|---|
| Product Code | CSB-EP013254BO |
| Storage | Store at -20°C, for extended storage, conserve at -20°C or -80°C. |
| Uniprot No. | P58754 |
| Product Type | Recombinant Protein |
| Immunogen Species | Bos taurus (Bovine) |
| Purity | >85% (SDS-PAGE) |
| Sequence | EPRE KMWICNSSDA SLWYNYCDDM KFPISVKVEP CVTIKGTKGK LHLYIARRD IQKLYLNLHI SIKSMTLPMR KEVICREYGG DYSFCGALKG ETVNTTIPFS FQGIRFSPGQ YHCVVEAISG NSEMLFCLN FTIIHYSSLN |
| Source | E.coli |
| Target Names | LY96 |
| Protein Names | Recommended name: Lymphocyte antigen 96 Short name= Ly-96 Alternative name(s): Protein MD-2 |
| Expression Region | 17-160 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Full Length of Mature Protein |
| Target Details | The MD-2 protein appears to associate with toll-like receptor 4 on the cell surface and confers responsiveness to lipopolysaccharide (LPS), thus providing a link between the receptor and LPS signaling. |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |