




Recombinant Rabbit Tropomyosin alpha-1 chain (TPM1)

| | |
|--------------------------|--|
| Product Code | CSB-YP024104RB |
| Relevance | Binds to actin filaments in muscle and non-muscle cells. Plays a central role, in association with the troponin complex, in the calcium dependent regulation of vertebrate striated muscle contraction. Smooth muscle contraction is regulated by interaction with caldesmon. In non-muscle cells is implicated in stabilizing cytoskeleton actin filaments. |
| Abbreviation | Recombinant Rabbit TPM1 protein |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | P58772 |
| Product Type | Recombinant Protein |
| Immunogen Species | Oryctolagus cuniculus (Rabbit) |
| Purity | ≥ 85% as determined by SDS-PAGE. |
| Sequence | MDAIKKKMQMLKLDKENALDRAEQAEADKKAEDRSKQLEDELVSLQKKLKG TEDELDKYSEALKDAQEKLELAEKKATDAEADVASLNRRIQLVEEELDRAQERL ATALQKLEEAEKAADESERGMKVIESRAQKDEEKMEIQEIQLKEAKHIAEDADR KYEEVARKLVIIESDLERAEEERAEELSEGKCAELEEEELKTVTNLKSLEAQAKEY SQKEDKYEEEIKVLSDKLKEAETRAEFAERSVTKLEKSIDDLEDELYAQKLYK AISEELDHALNDMTSI |
| Research Area | Signal Transduction |
| Source | Yeast |
| Target Names | TPM1 |
| Protein Names | Alpha-tropomyosin Tropomyosin-1 TPMA |
| Expression Region | 1-284aa |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | N-terminal 6xHis-tagged and C-terminal Myc-tagged |
| Mol. Weight | 36.2 kDa |
| Protein Length | Full Length |
| Image | |



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.