



# Recombinant Escherichia coli O17:K52:H18 ATP synthase subunit alpha (atpA)

<b>Product Code</b>	CSB-BP002344EOG
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	B7NF50
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli O17:K52:H18 (strain UMN026 / ExPEC)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MQLNSTEISE LIKQRIAQFN VVSEAHNEGT IVSVSDGVIR IHGLADCMQG EMISLPGNRY AIALNLERDS VGAVVMGPYA DLAEGMKVKC TGRILEVPVG RGLLGRVVNT LGAPIDGKGP LDHDGFSAVE AIAPGVIERQ SVDQPVQTGY KAVDSMIPIG RGQRELIIGD RQTGKTALAI DAIINQRDSG IKCIYVAIGQ KASTISNVVR KLEEHGALAN TIVVVATASE SAALQYLAPY AGCAMGEYFR DRGEDALIIY DDLISKQAVAY RQISLLLRP PGREAFPGDV FYLHSRLLER AARVNAEYVE AFTKGEVKGK TGSLTALPII ETQAGDVSAF VPTNVISITD GQIFLETNLF NAGIRPAVNP GISVSRVGG AQT KIMK KLS GGIRTALAQY RELAAFSQFA SDLDDATRQK LDHGQKVT EL LKQKQYAPMS VAQQSLV LFA AERGYLADVE LSKIGSFEAA LLAYVDRDHA PLMQEINQ TG GYNDEIEGKL KGILDSFKAT QSW
<b>Source</b>	Baculovirus
<b>Target Names</b>	atpA
<b>Protein Names</b>	Recommended name: ATP synthase subunit alpha EC= 3.6.3.14 Alternative name(s): ATP synthase F1 sector subunit alpha F-ATPase subunit alpha
<b>Expression Region</b>	1-513
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.