



# Recombinant Escherichia coli O7:K1 Porphobilinogen deaminase (hemC)

<b>Product Code</b>	CSB-YP010524EON
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	B7NTE0
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli O7:K1 (strain IAI39 / ExPEC)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MLDNVLR IAT RQSPLALWQA HYVKDKLMAS HPGLVVELVP MVTRGDVILD TPLAKVGGKG LFKVELEVAL LENRADI AVH SMKDVPEFEP QGLGLVTICE REDPRDAFVS NNYDNLDALP AGSIVGTSSL RRQCQLAERR PDLIIRSLRG NVGTRLSKLD NGEYDAIILA VAGLKRLGLE SRIRAALPPE ISLPAVGQGA VGIECRLDDT RTRELLAALN HHETALRVTA ERAMNTRLEG GCQVPIGSYA ELIDGEIWLR ALVGAPDGSQ IIRGERRGAP QDAEQMGISL AEELLNNGAR EILAEVYNGD APA
<b>Source</b>	Yeast
<b>Target Names</b>	hemC
<b>Protein Names</b>	Recommended name: Porphobilinogen deaminase Short name= PBG EC= 2.5.1.61 Alternative name(s): Hydroxymethylbilane synthase Short name= HMBS Pre-uroporphyrinogen synthase
<b>Expression Region</b>	1-313
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.