



Recombinant Escherichia fergusonii Peptide chain release factor 1 (prfA)

Product Code	CSB-MP015197EOR
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B7LSH9
Product Type	Recombinant Protein
Immunogen Species	Escherichia fergusonii (strain ATCC 35469 / DSM 13698 / CDC 0568-73)
Purity	>85% (SDS-PAGE)
Sequence	MKPSIVAKLE ALHERHEEVQ ALLGDAQTIA DQERFRALSR EYAQLSDVSR CFTDWQQVQE DIETAQMMLD DPEMREMAQD ELREAKEKSE QLEQQQLQVLL LPKDPDDERN AFLEVRAGTG GDEAALFAGD LFRMYSRYAE ARRWRVEIMS ASEGEHGGYK EIIAKISGDG VYGR LKFESG GHRVQRVPAT ESQGRIHTSA CTVAVMPELP DAELPDINPA DLRIDTFRSS GAGGQHVNTT DSAIRITHLP TGIVVECQDE RSQHKNKAKA LSVLGARIHA AEMAKRQQAE ASTRRNLLGS GDRSDRNRTY NFPQGRVTDH RINLTLYRLD EVMEGKLDML IEPIIQEHQA DQLAALSEQE
Source	Mammalian cell
Target Names	prfA
Protein Names	Recommended name: Peptide chain release factor 1 Short name= RF-1
Expression Region	1-360
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.