



# Recombinant Escherichia coli O17:K52:H18 3-ketoacyl-CoA thiolase (fadI)

<b>Product Code</b>	CSB-EP001117EOG-B
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	B7N5V3
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli O17:K52:H18 (strain UMN026 / ExPEC)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MGQVLPLVTR QGDRIAIVSG LRTPFARQAT AFHGIPAVDL GKMVVGELLA RSEIPAEVIE QLVFGQVVQM PEAPNIAREI VLGTGMNVHT DAYSVSRACA TSFQAIANVA ESLMAGTIRA GIAGGADSSS VLPIGVSKKL ARVLVDVNKA RTMSQRLKLF SRLRLRDLMP VPPAVAEYST GLRMGD TAEQ MAKTYGITRE QQDALAHRSH QRAAQAWSDG KLKEEVMTAF IPPYKQPLAE DNNIRGNSSL ADYAKLRPAF DRKHGTVTAA NSTPLTDGAA AVILMTESRA KELGLVPLGY LRSYAFTAID VWQDMLLGPA WSTPLALERA GLTMSDLTLI DMHEAFAAQT LANIQLLGSE RFARDVLGRA HATGEVDDSK FNVLGGSIA Y GHPFAATGAR MITQTLHEL RRGGGFGLVT ACAAGGLGAA MVVEAE
<b>Source</b>	E.coli
<b>Target Names</b>	fadI
<b>Protein Names</b>	Recommended name: 3-ketoacyl-CoA thiolase EC= 2.3.1.16 Alternative name(s): ACSs Acetyl-CoA acyltransferase Acyl-CoA ligase Beta-ketothiolase Fatty acid oxidation complex subunit beta
<b>Expression Region</b>	1-436
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.